

High-throughput Identification of Point Mutations in the Mouse Genome using Temperature Gradient Capillary Electrophoresis (TGCE)

Cymbeline T. Culiat, Ph. D.
Functional Genomics Group
Life Sciences Division


Functional Annotation of Human Genome Sequence

- Assigning function to ~30,000 - 50,000 genes
- Identifying regulatory elements (unknown number)
- **Experimentally** validating predictions of gene function
- Detailed analysis of mechanisms of function (complex biological pathways)

Mouse Mutagenesis: An Indispensable Approach

- Examine function by observing consequences of disrupting a gene in a **mammalian** system
- The mouse remains the premier mammalian animal model
- Mutant mouse serves as **reagent for both basic and pharmaceutical studies**
- Often yields **very surprising results!**

GABRB3 in Craniofacial Development



- GABRB3, a neurotransmitter receptor molecule was initially predicted to be a gene playing a major role in mental retardation
- The major role of this gene in facial clefting was first demonstrated at ORNL using a series of mouse mutations and gene-rescue experiments (Culiati *et al.*, 1995, Nature Genetics) .
- A recent publication (Scapoli *et al.*, 2002, Human Genetics) has **linked the GABRB3 gene to non-syndromic cleft lip and palate in man.**

Generating Mouse Mutations with ENU



William B Russell *et al.* (1979)
PNAS 76:5918-5911.
ORNL

- *N-ethyl-N-nitrosourea* (ENU), an alkylating agent
- Potent mutagen of mouse spermatogonial stem cells
- Induces primarily single base pair changes in DNA
- Different mutations in the same gene with a range of consequences- an allelic series
- Slightly altered proteins for protein-complex analysis

Large- scale ENU-mutagenesis Programs

- **Publicly funded**

- USA (ORNL, Jackson Lab, Baylor Univ, NIH consortia)
- United Kingdom (MRC, Harwell)
- Germany
- Japan
- Canada
- Australia

- **Privately funded**

- Novartis
- GlaxoSmithKline
- Pharmacia
- Ingenium Pharmaceuticals AG
- Hypnion Inc.
- Celltech R&D


Temperature-gradient capillary electrophoresis (TGCE)

Q. Gao and E. S. Yeung

High-throughput detection of unknown mutations by using multiplexed capillary electrophoresis with poly(vinylpyrrolidone) solution. *Analytical Chemistry* 2000, 72: 2499-2506

Ames National Laboratory, USDOE
Dept of Chemistry, Iowa State University

Mutation Scanning using the SCE9610 (SpectruMedix LLC)




- **High-throughput**
 - twelve 96-well plates/day
 - 150-800 bp PCR fragments with different Tms
- **Highly sensitive**
 - low sample requirement (~200 pg/ μ l)
 - 97% accuracy, 3% false positive rate
- **Simple and low-cost sample preparation**
 - mouse tail DNA preps for PCR templates
 - unpurified PCR reactions for TGCE scanning
 - ethidium bromide laser-induced fluorescence
- **Multiplexing**
 - different amplicons in single capillary
 - pooling of samples


Li, Q., Liu, Z., Monroe, H. and Culiati, C. T.
(2002) *Electrophoresis* 23:1499-1511

Mutation Scanning with TGCE

Heteroduplex formation




Separation by TGCE




- Retarded mobility of heteroduplexes compared to homoduplexes
- Optimal temperature for resolution is sequence-dependent


Optimizing TGCE




Fah 5961SB ; 493 bp, exon 7, G to A





Fah 6287SB ; 340 bp, exon 6, A to G



Regional Mutagenesis: Chromosome 7 p region



Craniosynostosis: Premature Fusion of Skull Bones



- Intracranial pressure
- Impaired cerebral flow
- Impaired vision and hearing
- Learning disabilities
- Psychological defects

A Gene-Driven Approach of Mouse Mutagenesis

Gene Selection

A set of genes in any biological pathway of interest

- Skin Carcinogenesis
- Craniofacial Birth Defects
- Low-dose Radiation Susceptibility
- Chemical and Pathogen Susceptibility
- Intracellular Transport
- Aging and DNA repair
- Maternal Effects on DNA Repair

Mutation Scanning

*CMMB: Cryopreserved Mutant Mouse Bank

A bank of DNAs, sperm and tissues from mice mutagenized by ENU.

Rederivation

Recover mutant mouse by IVF and/or ICSI


Functional Studies

Genetics
Pathology
Physiology
Biochemistry

*Initial Funding: LDRD

PIs: E. Michaud, C. Culiati and E. Rinchik

Generating the CMMB Resource




CMMB Progress

- First 2780 (55%) mice have been banked
- Established conditions for high-throughput PCR and TGCE mutation scanning
- **A 25-Mb screen under way; largest gene-driven mutation scanning experiment attempted to date.**
- Re-derivation from frozen sperm set-up
- FileMaker Pro Database

Mouse Mutagenesis at ORNL


**Regional
Phenotype-driven
Chrs 7, 10, 15, and X**

Aberrant Phenotype → Gene



**Whole-genome
Gene- and sequence-
driven**

Gene → Aberrant Phenotype



TGCE is a breakthrough in finding ENU-induced mutations in the mouse genome!

Acknowledgments

ORNL Staff

Ed Michaud
Bem Culiati
Dabney Johnson
Gene Rinchik
Mitch Klebig
Jay Snoddy
Frank Larimer

SpectruMedix LLC

Qingbo Li
Zhaowei Liu
Heidi Monroe
ZY Guo
Kevin Gutshall

ORNL Technical Staff

K. T. Cain	Sarah Shinpock
Debra Carpenter	Ann Wymore
Lori Easter	Marilyn Kerley
Alysyn Gardner	Lori Hughes
Kay Houser	Shen Lu
Bob Olszewski	Melissa York
Irina Pinn	Carmen Foster
Beverly Stanford	
Ginger Shaw	

DOE-BER



Marvin Frazier

LDRD Program



heterozygotes

controls



200bp, G>A @ 34bp





200bp, T>C @ 66bp





200bp, G>A @ 34bp &
T>C @ 66bp





493bp, G>A @ 268bp



341bp, G>C @ 98bp



341bp, G>C @ 98bp



A Crude Estimate of Mutation Frequency in the CMMB

- **1 mutation per 0.5 -2.5 Mb**
 - whole animal mutagenesis (Russell, Rinchik, ORNL)
 - mutagenized ES cell bank (T. Magnuson, UNC)
 - whole animal mutagenesis (S. Brown, Harwell)
- **An “average” gene (2.5 kb DNA) will have ~5 to 25 independent ENU-induced mutations**
- **Single SCE9610 can screen an amplicon/ week in entire bank**